

MEDICAL SEATING



Hausted Exam Furnishings offers a variety of medical stools to meet the needs of an active facility.

Manufactured with the highest quality materials to ensure a long life, Hausted stools offer the most wanted features in an exam stool.

- 16 in / 41 cm seat diameter to comfortably fit most users
- Pneumatic assist height adjustment for ease of use
- Choice of 3 base styles to complement facility décor
- Choice of height adjustment lever options
- Wide diameter base for greater stability
- Backrest and foot ring options
- 2 in / 5 cm Hard wheel casters
- 24 color choices
- DEHP-free vinyl coverings
- 350 lb / 159 kg weight capacity (EVENLY DISTRIBUTED)

9501 Series

Height: 18.5 in / 47 cm to 24 in / 61 cm

Base Diameter: 24 in / 61 cm

Seat Diameter: 16 in / 41 cm

Overall Height with Backrest: 33 in / 84 cm to 39 in / 99 cm

Weight Capacity: 350 lbs / 159 kg (EVENLY DISTRIBUTED)

Shipping Weight: 20 lbs / 9 kg



9501 Series is our most popular Physician exam stool – long lasting quality at an economical price.

9501-XX

Physician Stool With Single Lever Release and Black Composite Base

9501-B-XX

With Backrest

9501-FR-XX

With Foot Ring (18 in / 46 cm Diameter)

9501-BFR-XX

With Backrest and Foot Ring

9501-10BFR-XX

Lab/Counter Stool with Backrest and Foot Ring
Height Range: 24.5in / 62cm to 35in / 89cm

9702 Series

Height: 18.5 in / 47 cm to 24 in / 61 cm

Base Diameter: 24 in / 61 cm

Seat Diameter: 16 in / 41 cm

Overall Height with Backrest: 33 in / 84 cm to 39 in / 99 cm

Weight Capacity: 350 lbs / 159 kg (EVENLY DISTRIBUTED)

Shipping Weight: 23 lbs / 10 kg



9702-XX

Physician Stool with D-Handle Release and Polished Aluminum Base

9702-B-XX

With Backrest

9702-FR-XX

With Foot Ring (18 in / 46 cm Diameter)

9702-BFR-XX

With Backrest and Foot Ring

9702-10BFR-XX

Lab/Counter Stool with Backrest and Foot Ring
Height Range: 23.5in / 60cm to 34in / 86cm

9903 Series

Height: 18.5 in / 47 cm to 24 in / 61 cm

Base Diameter: 24 in / 61 cm

Seat Diameter: 16 in / 41 cm

Weight Capacity: 350 lbs / 159 kg (EVENLY DISTRIBUTED)

Shipping Weight: 23 lbs / 10 kg



9903-XX

Physician Stool with Chrome 360° Hand Release and Polished Aluminum Base

9903-FR

With Foot Ring (18 in / 46 cm Diameter)

9903-10FR-XX

Lab/Counter Stool with Foot Ring
22in / 56cm to 32in / 81cm

Standard Color Options

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Note: Due to the many variations in monitors and browsers, color samples may appear different. Swatch Card is available from your representative.



1.770.368.4700

Information contained herein is subject to change. The most current and complete product information can be found on our website.
www.hausted.com

Follow us on



© 2018, GF Health Products, Inc. All Rights Reserved. Hausted is a trademark of GF Health Products, Inc. is an ISO 13485 Certified Company.
GF1800080RevB18

graham field