



SOUNDCARE PLUS®

The SoundCare Plus is a professional ultrasound device, designed by clinicians. With this unique design, the SoundCare Plus comes equipped with two applicators, 1cm² and 5cm², that are plugged into the device allowing the user to switch between applicators with a single button push.

The SoundCare Plus has 20 presets and operates in 1 MHz and 3 MHz frequencies. The portable design lets you take it to your patient or keep it at the clinic. It's easy to use, and up and running within minutes.

PRODUCT FEATURES

- For the Treatment of Chronic and Acute Muscular Pain
- Includes Two Applicators, 1cm² and 5cm²
- Fully Functional 1 MHz and 3 MHz Frequencies
- Simultaneously Switch Between Applicators with the Push of a Single Button
- 20 User-Defined Presets with Pulsed and Continuous Therapy
- Variable Duty Cycle Outputs
- Ergonomic Handle Design
- Two-Year Warranty

ORDERING DETAILS

DQ9275 SoundCare Plus Professional Ultrasound Device

STANDARD ACCESSORIES

DQ9275W5 5 cm Applicator
 DQ9276MGC Medical Grade Power Cord
 LS5288 TheraSonic Ultrasound Gel, 8 oz
 DQ9275W1 1 cm² Applicator

TECHNICAL SPECIFICATIONS

Power 100 / 240 VAC
 Frequency 50 / 60 Hz
 US Frequency 1 & 3 MHz
 US Output Continuous & Pulsed