



SUPERSTIM PREMIUM ELECTRODES

The SuperStim premium electrode line is manufactured to have high tensile-strength lead wires. This is thanks to Leadstay, a proprietary bonding and moisture barrier layer protecting from delamination, particularly when used with moist heat therapy. A premium hydrogel provides comfort and uniform stimulation.

PRODUCT FEATURES:

- Available in Cloth or Foam and a Variety of Shapes and Sizes
- Strong Lead Wire Pull Strength
- Premium Dual Layer 35 mil Gel
- Gel Formulation Produces Maximum Dry-Out Resistance
- Optimum Distribution of Electrical Current
- Proudly Made in the USA

ORDERING DETAILS

Cloth

400-895	2" (5 cm) Round, 4/pk
400-899	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-893	3" (7.6 cm) Round, 4/pk
400-864	1.25" (3 cm) Round, 4/pk
400-858	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-889	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-854	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk
400-866	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

Foam

400-894	2" (5 cm) Round, 4/pk
400-898	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-892	3" (7.6 cm) Round, 4/pk
400-869	1.25" (3 cm) Round, 4/pk
400-857	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-888	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-853	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk
400-868	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk